

MRI - Protocol – Twin-CMR

Twin 1



Twin 2



Twin 3

3x preMOLL I 5(3)3 SAX Breath-hold	3x PSIR MOCO LGE SSFP SAX Breath--hold	(Stress) Prohance 0,1 mmol/Kg	Overview multi-slice bSSFP- PSIR.	post_MOLL I 4(1)3(1)2 SAX Breath-hold	PSIR MOCO LGE SSFP SAX Breath-hold	Prohance 0,05 mmol/Kg	Cine SAX- Paket LGE MOCO- PSIR SAX-stack 7/0 4CV, 3CV, 2CV	post_MOLL I 4(1)3(1)2 SAX Breath-hold	PSIR MOCO LGE SSFP SAX Breath-hold
--	---	--	--	--	--	-----------------------------	---	--	--

TI = 400ms

TI = 340

